



*An Efficient View to Life for More Than Seeing*

**ANKA 121-DM LED**  
*Mobile Operating Light*





### Head Specifications

Illumination (Lux)	120.000
Focusing Type	Electronic LED On/Off
Color Temperature Range	3.600-5.000 K
Endo Mode with Adjustable Intensity	Yes
CRI R <sub>a</sub>	> 96
CRI R <sub>g</sub>	> 97
Illumination Depth (L1+L2)	>203 cm
LED Life	60.000 hours
Radiant Energy Ratio	< 3,6 mW/m <sup>2</sup> Lux
Spot Size	23 - 31 cm

### System Specifications

Line Voltage	220 V 50 Hz $\pm$ % 10
Battery With Automatic Charge	Optional
Total Illumination (Lux)	120.000

### Accessories

Sterilizable Handle	Yes
---------------------	-----



Etkin Tıbbi Cihazlar Servis San. ve Ltd. Şti.

Address:1358 Sokak No: 9 / A Alsancak - İZMİR

Tel: (0 232) 465 25 41 - 464 16 57 - 464 00 20

Fax: (0 232) 464 14 93

etkintibbi@etkintibbi.com www.etkintibbi.com