

U0003

FUNCTIONS



Backrest adjustment : 0 - 68°



Thigh rest adjustment : 0 - 27°



Foot rest mechanical adjustment : 0 - 22°



Vascular position



Fowler position



Cardiac position



Trendelenburg position : 0 - 16°



Anti - trendelenburg position : 0 - 16°



Mattress platform minimum height : 49 cm (150 mm with castor)



Mattress platform maximum height : 79 cm (150 mm with castor)

OPTIONS

GROUP 1 / CASTOR OPTIONS



GROUP 3 / ELECTRICAL ACCESSORY OPTIONS



GROUP 4 / BACKREST MOTOR OPTIONS



GROUP 5 / LEGREST MOTOR OPTIONS











GROUP 6 / CONTROL PANEL OPTIONS



GROUP 7 / MATTRESS PLATFORM OPTIONS



GROUP 8 / OTHER ACCESSORY OPTIONS



DESCRIPTION OF OPTION CODES

- T16** CASTOR TWO WITH BRAKE (125 MM)
- T14** CASTOR CENTRAL & DIRECTIONAL LOCK (125 MM)
- T19** CASTOR CENTRAL & DIRECTIONAL LOCK (150 MM)
- K0** WITHOUT SIDE RAIL
- K5** WITH PE SIDE RAIL
- E0** WITHOUT ELECTRICAL ACCESSORIES
- E1** BATTERY RECHARGEABLE LEAD GEL TYPE. OUTPUT VOLTAGE: 24 V DC. CAPACITY: 1,2 Ah.
- E2** NURSE CONTROL PANEL. PROTECTION CATEGORY : IP 66. 15 FUNCTION KEY. UNDER BED LIGHTING SYSTEM
- E3** FOOT CONTROL PEDAL. ONE PEDAL CPR. ONE PEDAL EXAMINATION POSITION FROM EACH SIDE.
- E4** BACKREST MOTOR (4000 N)
- M1** BACKREST CPR MANUAL CONTROL MOTOR (2000 N)
- M2** LEGREST MOTOR (4000 N)
- M3** LEGREST CPR MANUAL CONTROL MOTOR (2000 N)
- M4** LEGREST MOTOR (4000 N)
- C1** HAND REMOTE CONTROL. PROTECTION CATEGORY : IP 66
- C2** SIDE RAIL CONTROL PANEL RIGHT
- C3** SIDE RAIL CONTROL PANEL LEFT
- Y1** REMOVABLE 4 PARTED MESH PLATFORM.
- Y3** REMOVABLE DOUBLE WALLED PE BED PLATFORM
- W0** WITHOUT OTHER ACCESSORIES
- W1** IV STAND : ONE HEIGHT ADJUSTMENT SYSTEM
- W2** LIFTING POLE : BEARING CAPACITY 75 Kg
- W3** X-RAY AVAILABLE FOR BACK SIDE WITH CASSETTE HOLDER AND COMPACT LAMINATE
- W4** TELEVISION STAND FOR 17" LCD TV
- W5** TRACTION APPARATUS
- W9** BED EXTENSION : 20 cm.
- W11** RESTRAINT STRAP

Must select only one of them.

Do not select the others if you select without accessory. If you select the others you can select all or 1-2-3 ... of them.

Must select minimum one. Can select 1-2-3 ... both.

</



* RL : Laminate color options page 82.