

i6

Productive .Powerful .Popular





CHISON i6 an exceptional cart-based color doppler system with sensitive touch panel, helps to reduce user fatigue, creates the indeed clinical value for physician and sonographers, who make contribution to deliver the best patient care and precision diagnose at the fingertips.

OPTIMIZED ERGONOMIC

-- Efficiently at your fingertips

Easy start from TOUCH Panel to finish day exams.

+ 19 inch LED monitor
10.4 inch touch panel

+ Keyboard
with backlit

+ CD / DVD

+ 4 Probe
Connectors



Touch panel

- Select preset
- Image adjustment
- Measure & Comment & Body mark
- Generate Report
- Image review



Wide Range Health




Powerful 4D

Innovation technologies, surprised imaging,
No limited light source, 360 degree.

Curved Panoramic

Curved Panoramic Imaging utilizes pattern recognition and image synthesis to generate wide and flexible view to reveal more anatomic information for diagnosis.




Quadplex


Quadplex combines B, color, and PW with automatic trace and measurement to help the user make diagnosis more conveniently and more accurately. Moreover, the three modes and measurement are updated in real time.

Elastography

- Display the elasticity of different tissues in different color.
- Provide more clinical information, especially for breast tumor, thyroid, liver and prostate.
- Strain ratio measurement quantitatively gives the ratio between the average strain of the selected region and of the nearby normal tissue region.




hcare Solution



Multi-follicle measurement

- Measure up to 16 follicles each time
- Quick and accurate measurement
- Convenient for the report management
- Improve the IVF successful rate



Continuous Wave (CW) Doppler

Detect blood flow with high velocity and help doctors diagnose with more clinical information.



4D/OB/GYN



ABD



Neonatal



MSK



Urology



Small parts



Vascular






General imaging

EXCEPTIONAL WORKFLOW

-- Focused on diagnosis

- Easyview™ for convenient review
- User-defined functions to fit your own work style
- Multiple output for quick data transfer (DICOM, DVI, VIDEO, VGA)
- Streamlined workflow for easy exam completion







Ovary, 2B Mode



Fetal Face, Virtual HD




Hepatic Vein, PIP Zoom




Fetus, B/BC Mode




Early Pregnancy, C Mode



Lymph Node Enlargement, B Mode



Testicle, CPA Mode



Carotid Vessel, C Mode



2.5MHz-5.3MHz Convex
D3C60L



5.0MHz-10.0MHz Linear
D7L40L



4.0MHz-10.0MHz Transvaginal
D6C12L



4.0MHz-10.7MHz Transvaginal
D7C10L



2.5MHz-5.3MHz Volume
V4C40L



2.5MHz-6.4MHz Micro-Convex
D3C20L



2.5MHz-4.0MHz Phased array
D3P64L



4.0MHz-10.7MHz Micro-Convex
D6C15L



3.5MHz-8.0MHz Micro-Convex
D5C20L

CHISON Medical Technologies Co., Ltd.

Sales & Service Contact Address:

No.9, Xinhuihuan Road, Xinwu District, Wuxi, Jiangsu, China 214028

TEL : 0086-510-85310593 / 85310937

FAX : 0086-510-85310726

EMAIL : export@chison.com.cn

We reserve the right to make changes to this catalogue without prior notice
Please contact our local dealer for the latest information.

i6-20171113