

# STABILUS



## Product Feature

Features:

- \*Compact and durable. Adaptable to various working environment. Cost effective with high quality.
- \* High-quality medical-class stainless steel. Specially treated matte surface. Antibacterial and easy cleaning.
- \*Integrated waist bridge can be lifted up to 110 mm , which meets the needs of urinary surgery and other occasions.
- \*Newly designed self-lubricant lifting column. No need to fill lubricant oil.
- \*Heavy-duty casters with low noise and high mobility.
- There are standard table top and sliding table top available.
- \*“Screw-down” plinth stopper for a safe and steady operation.

\*Reliable hydraulic system for a steady and swift lifting up/down.

\*Pneumatic table adjustment system. Very easy and light to set to any position.

\*High quality PU surface pad. Ultrasound welding, no stitch leakage . Antistatic and chemical environment durable.

\*Dual-layer slow-elastic memory foam pad. More comfort and balanced stress than normal memory foam pad.

\*Orthopedic Tractor and Skull clamp as optional.

\*SFDA and CE approved

Pistons and Casters made in Germany

Assembled in Egypt